

Mass Spectrum List Report

Analysis Info


Analysis Name D:\Data\ggc\LR-MM-5_P1-C-1_01_8334.d
Method lc apci directos.m
Sample Name LR-MM-5
Comment Sv: DCM + MeOH + HCOOH
Durán

Acquisition Date 12/6/2023 10:31:38 AM

Operator ELB
Instrument / Ser# microTOF-Q II

Acquisition Parameter

Source Type	APCI	Ion Polarity	Positive	Set Nebulizer	3.0 Bar
Focus	Not active	Set Capillary	4500 V	Set Dry Heater	200 °C
Scan Begin	100 m/z	Set End Plate Offset	-500 V	Set Dry Gas	6.0 l/min
Scan End	1300 m/z	Set Collision Cell RF	250.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	237.0758	49341	22.8
2	238.0761	10596	4.9
3	281.0453	10814	5.0
4	299.0571	11562	5.3
5	311.0928	11011	5.1
6	313.2683	40954	18.9
7	314.2709	9200	4.3
8	329.1696	9139	4.2
9	341.2969	30452	14.1
10	355.0639	52319	24.2
11	356.0640	17285	8.0
12	357.0619	11313	5.2
13	371.0939	22356	10.3
14	372.0949	8752	4.0
15	429.0817	45331	21.0
16	430.0803	18638	8.6
17	431.0797	13365	6.2
18	445.1115	14954	6.9
19	503.0990	34019	15.7
20	504.0984	15603	7.2
21	505.0973	12250	5.7
22	519.1313	15917	7.4
23	551.4955	180074	83.3
24	552.4982	67320	31.1
25	553.4994	12541	5.8
26	579.5263	216228	100.0
27	580.5291	82602	38.2
28	581.5308	17042	7.9
29	607.5562	59409	27.5
30	608.5581	23542	10.9