

Mass Spectrum List Report

Analysis Info

Analysis Name D:\Data\ggc\24-110803_P1-E-2_01_9721.d
 Method lc-qtof directos tune high 141124 pos.m
 Sample Name 24-110803
 Comment Muestra 2
 Sv: MeOH + Agua
 Iribarren


Acquisition Date 11/14/2024 3:23:36 PM

Operator EB - MV

Instrument / Ser# microTOF-Q

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	2.0 Bar
Focus	Not active	Set Capillary	3000 V	Set Dry Heater	200 °C
Scan Begin	400 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	3100 m/z	Set Collision Cell RF	700.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	531.3500	14253	44.8
2	575.3780	22602	71.1
3	591.3527	12722	40.0
4	619.4029	22837	71.8
5	620.4061	7257	22.8
6	635.3810	24172	76.0
7	636.3851	7904	24.9
8	663.4318	20810	65.4
9	679.4080	28615	90.0
10	680.4094	9164	28.8
11	707.4540	13124	41.3
12	723.4322	23528	74.0
13	724.4358	8197	25.8
14	729.4969	6960	21.9
15	751.4807	7490	23.6
16	757.5430	9789	30.8
17	767.4543	15332	48.2
18	773.5195	18755	59.0
19	774.5250	7722	24.3
20	801.5701	12290	38.6
21	811.4831	8921	28.1
22	817.5462	29885	94.0
23	818.5507	12308	38.7
24	845.5957	13243	41.6
25	861.5759	31802	100.0
26	862.5723	13534	42.6
27	889.6228	7931	24.9
28	905.6008	24457	76.9
29	906.5983	11294	35.5
30	949.6299	13780	43.3