

Mass Spectrum List Report

Analysis Info

Analysis Name D:\Data\ggc\2Escuarato 30-9-24_P1-B-1_01_10429.d
 Method lc_ms directos tunelow291025.m
 Sample Name 2Escuarato 30-9-24
 Comment Sv: DCM + MeOH + HCOOH
 Marino


Acquisition Date 11/4/2025 2:24:28 PM

Operator EB - MV

Instrument / Ser# microTOF-Q

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	2.0 Bar
Focus	Not active	Set Capillary	4000 V	Set Dry Heater	200 °C
Scan Begin	100 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	1500 m/z	Set Collision Cell RF	250.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	149.0350	1031	7.8
2	225.0529	1507	11.3
3	226.9618	3268	24.6
4	233.1756	3652	27.5
5	259.1644	943	7.1
6	273.1777	2999	22.6
7	281.0629	1913	14.4
8	299.0770	4362	32.8
9	300.0769	1196	9.0
10	301.1561	2326	17.5
11	304.3108	863	6.5
12	315.1947	1045	7.9
13	325.2430	3941	29.7
14	326.2462	1008	7.6
15	331.2189	1442	10.9
16	341.2815	1280	9.6
17	353.2768	2244	16.9
18	355.0826	2483	18.7
19	381.3142	2220	16.7
20	385.3041	1072	8.1
21	393.3114	2144	16.1
22	413.2832	13284	100.0
23	414.2879	3168	23.8
24	415.2924	1180	8.9
25	427.2654	1894	14.3
26	441.3173	4899	36.9
27	442.3196	1535	11.6
28	503.1281	1455	11.0
29	553.4819	2711	20.4
30	610.2088	1393	10.5