

Mass Spectrum List Report

Analysis Info


Analysis Name D:\Data\ggc\HG-161_P1-D-2_01_9226.d
 Method lc_qtof_directospos.m
 Sample Name HG-161
 Comment Sv: MeOH + HCOOH
 Szajnman

Acquisition Date 6/6/2024 10:52:08 AM

Operator MV
 Instrument / Ser# microTOF-Q II

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	3.5 Bar
Focus	Not active	Set Capillary	4000 V	Set Dry Heater	200 °C
Scan Begin	70 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	1000 m/z	Set Collision Cell RF	300.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	197.0603	16615	15.0
2	219.0412	111102	100.0
3	220.0436	13771	12.4
4	226.9467	6551	5.9
5	229.1362	9697	8.7
6	266.0658	7991	7.2
7	281.0535	26167	23.6
8	283.0685	10055	9.1
9	299.0660	9793	8.8
10	301.1320	27351	24.6
11	329.0708	7416	6.7
12	353.2564	13726	12.4
13	358.3587	11401	10.3
14	365.1295	8277	7.4
15	367.1249	15663	14.1
16	369.1341	6981	6.3
17	379.2686	6741	6.1
18	381.2901	12992	11.7
19	413.2567	17199	15.5
20	509.1418	70203	63.2
21	510.1475	20185	18.2
22	511.1516	7562	6.8
23	525.1361	15191	13.7
24	527.1510	7090	6.4
25	553.4443	9201	8.3
26	661.2091	6651	6.0
27	685.4205	14166	12.8
28	807.1437	7935	7.1
29	809.1473	14987	13.5
30	810.1476	7377	6.6