

Mass Spectrum List Report

Analysis Info


Analysis Name D:\Data\ggc\HG-163_P1-D-4_01_9228.d
 Method lc_qtof_directospos.m
 Sample Name HG-163
 Comment Sv: MeOH + HCOOH
 Szajnman

Acquisition Date 6/6/2024 11:02:22 AM

Operator MV
 Instrument / Ser# microTOF-Q II

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	3.5 Bar
Focus	Not active	Set Capillary	4000 V	Set Dry Heater	200 °C
Scan Begin	70 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	1000 m/z	Set Collision Cell RF	300.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	226.9515	8980	36.1
2	229.1413	8767	35.2
3	270.9787	4216	16.9
4	301.1416	11034	44.3
5	353.2632	11101	44.6
6	354.2670	3106	12.5
7	358.3667	15846	63.7
8	359.3616	4150	16.7
9	365.1244	3166	12.7
10	379.2772	5353	21.5
11	381.2961	11174	44.9
12	392.9758	5040	20.3
13	393.2921	3465	13.9
14	413.2648	24887	100.0
15	414.2662	6493	26.1
16	441.2897	2628	10.6
17	509.1509	2869	11.5
18	525.0306	4198	16.9
19	591.4925	5458	21.9
20	596.5975	3955	15.9
21	619.5232	6939	27.9
22	620.5256	3145	12.6
23	624.6211	5237	21.0
24	633.0333	3290	13.2
25	635.0499	4193	16.8
26	663.4569	5219	21.0
27	679.4287	2647	10.6
28	680.4616	2789	11.2
29	685.4286	5065	20.4
30	881.7459	3020	12.1