

Mass Spectrum List Report

Analysis Info

Analysis Name D:\Data\ggc\HG-3i_P1-C-6_01_9887.d
 Method lc-ms directos tunelow261224 apci.m
 Sample Name HG-3i
 Comment Sv: MeOH + HCOOH
 Szajnman


Acquisition Date 12/26/2024 4:12:11 PM

Operator EB - MV

Instrument / Ser# microTOF-Q

Acquisition Parameter

Source Type	APCI	Ion Polarity	Positive	Set Nebulizer	1.6 Bar
Focus	Not active	Set Capillary	4500 V	Set Dry Heater	200 °C
Scan Begin	100 m/z	Set End Plate Offset	-500 V	Set Dry Gas	8.0 l/min
Scan End	1000 m/z	Set Collision Cell RF	250.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	145.0567	8135	4.7
2	159.0979	5293	3.1
3	171.1031	16313	9.4
4	173.1173	11912	6.9
5	187.1352	84915	49.0
6	188.1374	9083	5.2
7	189.1336	4234	2.4
8	201.1474	5663	3.3
9	213.1474	10047	5.8
10	215.1660	173363	100.0
11	216.1687	20657	11.9
12	217.1624	7313	4.2
13	227.1639	20674	11.9
14	231.1585	48110	27.8
15	232.1619	6781	3.9
16	239.1588	3624	2.1
17	263.1596	19763	11.4
18	313.2663	77748	44.8
19	314.2718	13346	7.7
20	315.2683	3730	2.2
21	325.1935	14439	8.3
22	331.1901	4895	2.8
23	349.2253	3451	2.0
24	359.2485	7314	4.2
25	385.2643	8903	5.1
26	387.2740	7520	4.3
27	557.3746	19998	11.5
28	558.3853	5695	3.3
29	559.3757	3668	2.1
30	571.3951	3761	2.2