

# Mass Spectrum List Report

## Analysis Info


Analysis Name D:\Data\ggc\HG-3i\_P1-C-3\_01\_9988.d  
Method lc\_ms directos tunelow300924.m  
Sample Name HG-3i  
Comment Sv: MeOH  
Szajnman

Acquisition Date 1/15/2025 3:35:34 PM

Operator EB - MV  
Instrument / Ser# microTOF-Q

## Acquisition Parameter

|             |            |                       |           |                  |           |
|-------------|------------|-----------------------|-----------|------------------|-----------|
| Source Type | ESI        | Ion Polarity          | Positive  | Set Nebulizer    | 2.0 Bar   |
| Focus       | Not active | Set Capillary         | 4000 V    | Set Dry Heater   | 200 °C    |
| Scan Begin  | 100 m/z    | Set End Plate Offset  | -500 V    | Set Dry Gas      | 7.0 l/min |
| Scan End    | 1200 m/z   | Set Collision Cell RF | 250.0 Vpp | Set Divert Valve | Source    |



| #  | m/z      | I      | I %   |
|----|----------|--------|-------|
| 1  | 171.1134 | 2070   | 1.2   |
| 2  | 187.1477 | 10694  | 6.3   |
| 3  | 215.1788 | 37763  | 22.4  |
| 4  | 216.1806 | 4592   | 2.7   |
| 5  | 217.1664 | 3002   | 1.8   |
| 6  | 231.1739 | 168889 | 100.0 |
| 7  | 232.1748 | 20445  | 12.1  |
| 8  | 233.1716 | 7481   | 4.4   |
| 9  | 313.2862 | 6456   | 3.8   |
| 10 | 325.2140 | 63578  | 37.6  |
| 11 | 326.2166 | 11496  | 6.8   |
| 12 | 327.2128 | 3620   | 2.1   |
| 13 | 347.1954 | 45498  | 26.9  |
| 14 | 348.1975 | 7401   | 4.4   |
| 15 | 349.2028 | 2703   | 1.6   |
| 16 | 359.2711 | 2484   | 1.5   |
| 17 | 387.2986 | 2629   | 1.6   |
| 18 | 413.2590 | 3303   | 2.0   |
| 19 | 415.2932 | 2889   | 1.7   |
| 20 | 469.3188 | 3544   | 2.1   |
| 21 | 471.3255 | 3333   | 2.0   |
| 22 | 501.3314 | 9726   | 5.8   |
| 23 | 502.3358 | 2129   | 1.3   |
| 24 | 503.3320 | 2394   | 1.4   |
| 25 | 517.3465 | 2303   | 1.4   |
| 26 | 527.4020 | 2048   | 1.2   |
| 27 | 557.4096 | 4967   | 2.9   |
| 28 | 619.4083 | 2340   | 1.4   |
| 29 | 643.4422 | 2706   | 1.6   |
| 30 | 663.4320 | 2254   | 1.3   |