

Mass Spectrum List Report

Analysis Info


Analysis Name D:\Data\ggc\HG-3i_P1-C-6_01_9828.d
 Method lc_ms directos tunelow300924.m
 Sample Name HG-3i
 Comment Sv: MeOH + HCOOH
 Szajnman

Acquisition Date 12/23/2024 1:11:31 PM

Operator EB - MV
 Instrument / Ser# microTOF-Q

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	2.0 Bar
Focus	Not active	Set Capillary	4000 V	Set Dry Heater	200 °C
Scan Begin	100 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	1000 m/z	Set Collision Cell RF	250.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	187.1621	7132	14.9
2	215.1975	18429	38.4
3	231.1934	47964	100.0
4	232.1965	5504	11.5
5	279.1108	3179	6.6
6	313.3128	6604	13.8
7	325.2412	16629	34.7
8	326.2441	2940	6.1
9	344.2181	19070	39.8
10	344.7207	7662	16.0
11	345.2196	5210	10.9
12	347.2245	37884	79.0
13	348.2290	5688	11.9
14	353.2207	2922	6.1
15	355.0958	4113	8.6
16	365.2101	2790	5.8
17	393.2045	3553	7.4
18	413.2901	8124	16.9
19	429.1185	4394	9.2
20	469.3561	8003	16.7
21	471.3635	6312	13.2
22	501.3734	19040	39.7
23	502.3811	4514	9.4
24	503.3716	3541	7.4
25	506.3389	7355	15.3
26	506.8377	4694	9.8
27	517.3866	3015	6.3
28	557.4654	5965	12.4
29	619.4585	3372	7.0
30	643.4998	3413	7.1