

Mass Spectrum List Report

Analysis Info


Analysis Name D:\Data\ggc\V80p196_P1-B-6_01_9017.d
 Method lc_qtof_directospos.m
 Sample Name V80p196
 Comment Sv: MeOH + HCOOH
 Uhrig

Acquisition Date 4/18/2024 3:49:58 PM

Operator ELB
 Instrument / Ser# microTOF-Q II

Acquisition Parameter

Source Type	ESI	Ion Polarity	Positive	Set Nebulizer	3.5 Bar
Focus	Not active	Set Capillary	4000 V	Set Dry Heater	200 °C
Scan Begin	70 m/z	Set End Plate Offset	-500 V	Set Dry Gas	7.0 l/min
Scan End	1000 m/z	Set Collision Cell RF	300.0 Vpp	Set Divert Valve	Source



#	m/z	I	I %
1	192.0952	10756	4.3
2	206.0830	15542	6.2
3	208.0924	18294	7.2
4	214.1316	4831	1.9
5	222.1093	8410	3.3
6	226.9646	7365	2.9
7	229.0949	27899	11.1
8	234.1095	6645	2.6
9	261.1334	9591	3.8
10	275.1136	252380	100.0
11	276.1156	27608	10.9
12	277.1281	9587	3.8
13	316.1437	5309	2.1
14	318.1586	11307	4.5
15	332.1684	5065	2.0
16	334.1581	4587	1.8
17	337.2492	4478	1.8
18	350.1859	8303	3.3
19	360.3429	9127	3.6
20	361.1816	16444	6.5
21	362.1861	4549	1.8
22	373.1918	4948	2.0
23	375.1939	5417	2.1
24	387.1926	8592	3.4
25	389.1962	4470	1.8
26	425.1584	6298	2.5
27	458.2052	7076	2.8
28	481.2168	9096	3.6
29	527.2345	67335	26.7
30	528.2396	14960	5.9